

Protocol	#9.1
Title	BOMB bisulfite conversion
Keywords	DNA methylation, HT bisulfite conversion, Silica-beads
Authors	Oberacker P*, Stepper P*, Bond DM*, Höhn S, Focken J, Meyer V, Schelle L, Sugrue VJ, Jeunen
	GJ, Moser T, Hore SR, von Meyenn F, Hipp K, Hore TA# and Jurkowski TP#
Citation	Oberacker et al. (2019), Bio-On-Magnetic-Beads (BOMB): Open platform for high-throughput
	nucleic acid manipulation. PLOS Biology, 17(1), https://doi.org/10.1371/journal.pbio.3000107
Online	https://bomb.bio/protocols/
Revision	V1.0 (13 <sup>th</sup> August 2018)

## Summary

Bisulfite sequencing [1] is considered the gold standard [2] for the analysis of DNA methylation status. Purified genomic DNA is treated with sodium metabisulfite which converts unmodified cytosines to uracils (which later amplify as thymines). Methylated cytosines (5mC) are resilient to bisulfite conversion and amplify as Cs. Thereby, the methylation status of each CpG site can be determined by sequencing and comparison to an unconverted reference sequence [3,4]. This is a modified, fast bisulfite conversion protocol and includes a rapid conversion, magnetic separation of treated DNA and on-bead desulfonation. This protocol was optimized for DNA amounts of 500-800 ng, but works well with input amounts ranging from 10-20 ng to 2  $\mu$ g.

## Chemicals

Name	Provider	PN	MW [g/mol]	Safety codes	
Ethanol (C₂H <sub>6</sub> O, 99.9%)	Honeywell / Riedel-de Haën	34963	46.07	<b>anger</b>	H: 225-319 P: 210-80- 305+351+338- 308+313
Sodium metabisulfite (Na₂S₂O₅)	Sigma-Aldrich (Merck)	S9000-500G	190.11	<b>Danger</b>	H: 302-318 P: 280-301+312+330- 305+351+338+310
Hydroquinone (C₀H₀O₂)	Sigma-Aldrich (Merck)	H9003-100G	110.11	L Danger	H: 302-317-318-341- 351-410 P: 201-273-280- 305+351+338+310- 308+313
N,N- Dimethylformamide (DMF, C₃H7NO)	Acros Organics (Thermo Fisher Scientific)	210580025	73.09	<ul><li>Danger</li></ul>	H: 226-312-319-332- 360D P: 301+310-302+352- 280-304+340- 305+351+338-210
Guanidine hydrochloride (CH₅N₃ · HCl)	Roth Chemicals	0037.1	95.53	() Warning	H: 302-315-319 P: 280-302+352- 305+351+338
Sodium hydroxide (NaOH)	Roth Chemicals	6771.1	40.00	Danger	H: 290-314 P: 280-301+330+331- 305+351+338-310

Please consult appropriate MSDS information before working with these chemicals! Use lab coat, gloves and eye protection at all times! The chemicals are available from other providers as well. No preference is given to the indicated vendors.





# **Buffers and solutions**

Silica-coated magnetic beads (20  $\mu$ l per sample, 2 ml for a 96-well plate)

Take 1 ml of resuspended silica-coated magnetic beads (see **BOMB protocol #2.1**), pellet magnetically, and discard supernatant. Wash once with 1 ml of 7 M guanidine-HCl and resuspend in 2 ml 7 M guanidine-HCl

**Conversion reagent** (130  $\mu$ l per sample, 13 ml for a 96-well plate, prepare freshly)

Reagent	Per 96-well plate
Sodium metabisulfite (Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> )	6 g
Hydroquinone	5 mg
DMF	250 μl of stock
NaOH	3 ml of 2 M stock
ddH <sub>2</sub> O	9.25 ml

Desulfonation solution (prepared fresh. 20 ml per plate)

0.25 M NaOH (2.5 ml of 2 M) 80% ethanol (16 ml of absolute, 1.5 ml water)

Binding buffer (if not dissolving, heat to 35 °C. 60 ml per plate)

7 M guanidine-HCl

90% ethanol (250 µl per sample, 25 ml per 96) - can be stored at RT for at least 12 months if closed properly

**Elution buffer** (60 µl per sample, 6 ml per 96) – can be stored at RT for at least 12 months

5 mM Tris-HCl pH 8.5

## **Equipment and setup**

Fume hood

Heating block (e.g. Roth, Rotilabo®-block thermostat H250 – PN: Y264.1)

**Thermocycler** (e.g. BIOMETRA T1)

Microtiter-plate orbital shaker (e.g. IKA MS 3 basic)

Magnetic stand for 96-well plate (e.g. BOMB microplate magnetic rack)

#### **Multichannel Pipettes**

30 – 300 μl (e.g. Eppendorf, Eppendorf Research<sup>®</sup> plus 8-Channel – PN: 3125000052)

50 – 1200 µl (e.g. VWR, Multi-channel pipette, 8-channel – PN: 613-5422)

96-well PCR plate (e.g. Sarstedt - PN:72.1979.102)

1.2 ml 96-well deep well plates (e.g. Sarstedt, MegaBlock® 96 Well – PN: 82.1971.002)

Reservoirs (e.g. Roth, Rotilabo®-liquid reservoirs – PN: E830.2)

Seals (e.g. Bio-Rad, Microseal<sup>®</sup> 'B' Adhesive Seals – PN: MSB1001)



This work is licensed under a <u>Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License</u>



# **BOMB Bisulfite conversion**

Step	Bisulfite conversion	Time	$\mathbf{\overline{\mathbf{A}}}$
1	Transfer 20 ng to 1 μg of your DNA sample in a total volume of 20 μl to fresh 0.2 ml PCR tubes or 96-well PCR plate		
2	To each sample add 130 $\mu l$ of <b>conversion reagent</b> , seal, mix well and spin down	E min	
$\bigwedge$	DMF is highly toxic; sodium hydroxide and bisulfite are corrosive and dangerous!	5 11111	
3	Place the samples in a thermocycler with a heated lid and execute the following program: Step 0 – 95 °C pause Step 1 – 95 °C 5 min Step 2 – 54 °C 30 min Step 3 – 95 °C 1 min Step 4 – 54 °C 30 min Step 5 – 95 °C 1 min Step 6 – 54 °C 30 min Step 7 – 4 °C storage for up to 20 hours Converted DNA without further purification should not be stored for extended period of time due to progressing DNA degradation resulting in loss of DNA and decreased efficiency of the subsequent amplification step	2 h	
Step	Desulfonation and DNA purification		
4	Pipette 20 $\mu$ l of <b>silica-coated magnetic beads</b> (in 7 M guanidine-HCl) into the wells of a deep-well plate and add 600 $\mu$ l of <b>binding buffer</b>	5 min	
5	Transfer the <b>converted sample</b> (150 $\mu$ l from step 3) to the deep-well plate with beads in binding buffer	5 min	
6 ⁄!	Shake for 5 min to mix and bind If necessary, seal the plate and invert a few times to resuspend the beads	5 min	
7	Pellet the beads with a magnetic rack and remove the supernatant. Wash twice with 400 $\mu l$ of <b>90% ethanol</b>	10 min	
8	Desulfonate the DNA by adding 200 $\mu I$ of $desulfonation buffer$ and shaking at RT for $15\ min$	20 min	
9	Pellet the beads with a magnetic rack, remove the supernatant and wash with 300 $\mu l$ of <b>90% ethanol</b>	5 min	
10	Wash twice with 200 μl of <b>90% ethanol</b>	10 min	
11	Dry the beads by placing the plate at 50 °C for 10-30 min (until the beads become brittle and brown)	30 min	
12	Add 60 $\mu$ l of <b>elution buffer</b> to the wells and shake shortly to resuspend (centrifuge shortly if beads stick to the walls)	5 min	
13	Pellet the beads magnetically and pipette off the eluted <b>bisulfite converted DNA</b> into a fresh plate or tubes	5 min	
Fnd	Measure concentration	re concentration ~3.5 h	
		(1.5 h hai	nds-on)

Store @ -80 °C (several weeks @ 4 °C is preferable to freeze-thaw cycles)



# Troubleshooting

Problem	Solution
Beads sticking to the sides	• Sonicate shortly in a sonic bath and/or push the beads down with a pipette tip
Beads stay at the bottom of the well when mixing	<ul> <li>Use a good foil to seal the plate and invert a few times to mix. Spin down briefly</li> </ul>
Low conversion efficiency	Ensure that you're using fresh conversion solutions
Low DNA amplification	<ul> <li>DNA degradation: Use more input DNA. If that is not possible, test skipping one cycle of 95 and 54 °C during conversion in the thermocycler</li> <li>Suboptimal primers: use a set of validated PCR primers that are efficient in</li> </ul>

amplification of converted DNA

### **Exemplary Results**



**Figure S5: Bisulfite conversion.** (A) Scheme of DNA methylation analysis using bisulfite conversion. (B) Agarose gel after PCR amplification of bisulfite converted human DNA. Multiple primer pairs with expected product sizes between 221 and 435 bp were tested successfully. MW: GeneRuler DNA Ladder Mix, (Thermo Scientific). (C) Sequencing trace of a PCR amplified bisulfite converted sample aligned with the original, unconverted sequence. All non-CG cytosines were successfully converted. Conversion rate is ~99% as measured by sequencing after PCR amplification

# References

- 1. Frommer M, McDonald LE, Millar DS, Collis CM, Watt F, Grigg GW, et al. A genomic sequencing protocol that yields a positive display of 5-methylcytosine residues in individual DNA strands. Proc Natl Acad Sci U S A. 1992;89: 1827–31. doi:https://doi.org/10.1073/pnas.89.5.1827
- 2. Kurdyukov S, Bullock M. DNA Methylation Analysis: Choosing the Right Method. Biology (Basel). Multidisciplinary Digital Publishing Institute; 2016;5: 3. doi:10.3390/biology5010003
- Shiraishi M, Hayatsu H. High-Speed conversion of cytosine to uracil in bisulfite genomic sequencing analysis of DNA methylation. DNA Res. Oxford University Press; 2004;11: 409–415. doi:10.1093/dnares/11.6.409
- 4. Wang H, Ke H, Zheng Y, Lai J, Luo Q, Chen Q. A modified bisulfite conversion method for the detection of DNA methylation. Epigenomics. 2017;9: 955–969. doi:10.2217/epi-2016-0174