









| | |
|-----------------|---|
| Protocol | #5.4 |
| Title | BOMB plasmid MIDI prep using Sera-Mag carboxylated beads |
| Keywords | midiprep, Sera-Mag, plasmid extraction |
| Authors | Treitli S. C. |
| Citation | <i>Treitli S.C., BOMB.bio, 2019</i> |
| Online | https://bomb.bio/protocols/ |
| Editing | V1.0 (Oberacker P., 16 th September 2019) |

Chemicals

| Name | Provider | PN | MW [g/mol] | Safety codes | |
|---|----------------------------|----------|---------------|--|---|
| Acetic acid (CH₃COOH) | Merck | 1.09951 | 60.05 |  Danger | H: 290-314 P: 280-301+330+331-305+351+338-308+310 |
| Ethanol (C₂H₆O, 99.9%) | Honeywell / Riedel-de Haën | 34963 | 46.07 |  Danger | H: 225-319 P: 210-280-305+351+338-308+313 |
| Ethylenediaminetetraacetic acid dihydrate / EDTA (C₁₀H₁₆N₂O₈ · 2 H₂O) | Roth Chemicals | 8043.1 | 372.24 |  Warning | H: 332-373 P: 260-314 |
| Hydrochlorid acid fuming (HCl_(aq), 37%) | Roth Chemicals | 4625.2 | 36.46 |  Danger | H: 290-314-335 P: 280-303+361+353-304+340-305+351+338-310 |
| Polyethylene glycol 8000 | AppliChem | A2204 | ~8000 | n.a. | n.a. |
| Potassium acetate (CH₃CO₂K) | Riedel-de Haën | 32309 | 98.15 | n.a. | n.a. |
| RNase A | Serva | 34390.02 | n.a. |  Danger | H: 334 P: 261-284-304+340-342+311 |
| Sodium Chloride (NaCl) | Roth Chemicals | 3957.2 | 58.44 | n.a. | n.a. |
| Sodium dodecyl sulfate / SDS (C₁₂H₂₅NaO₄S) | Roth Chemicals | CN30.3 | 288.38 |  Danger | H 302-315-318-412 P: 280-301+312-302+352-305+351+338-332+313 |

| Name | Provider | PN | MW [g/mol] | Safety codes | |
|--|-------------------|--------|---------------|--|--|
| Sodium hydroxide (NaOH) | Roth Chemicals | 6771.2 | 40.00 |  Danger | H: 290-314 P: 280- 301+330+331- 305+351+338- 310 |
| Tris(hydroxymethyl)- aminomethane / Tris (C ₄ H ₁₁ NO ₃) | Roth Chemicals | AE15.3 | 121.14 |  Warning | H: 315-319-335 P: 280-302+352- 305+351+338- 312 |
| Tween 20 | Roth Chemicals | 9127.1 | 1228 | n.a. | n.a. |

Please consult appropriate MSDS information before working with these chemicals! Use lab coat, gloves and eye protection at all times! The chemicals are available from other providers as well. No preference is given to the indicated vendors.

Buffers and solutions

1 M Tris (50 ml)

6.057 g – TRIS

40 ml of MilliQ water

Adjust the pH to 8.0 with concentrated HCl

Bring the final volume to 50 mL

Filter sterilize.

0.5 M EDTA (50 ml)

9.31g – EDTA

40 ml of MilliQ water

Adjust the pH to 8.0 NaOH. You will need around one gram of NaOH. When you get close to pH 8, it is better to use NaOH solution, to avoid overshooting the pH.

Bring the final volume to 50 mL

Filter sterilize.

1 M NaOH (50 ml)

2g -NaOH

Dissolve in 40 ml of MilliQ water.

Bring the final volume to 50 ml

P1 buffer – can be stored at 4 °C for up to 12 months

| Reagent | Concentrations | For 50 ml |
|-----------------|----------------|-----------------------------|
| Tris-HCl pH 8.0 | 50 mM | 2.5 ml of 1 M stock |
| EDTA | 10 mM | 1 ml of 0.5 M stock |
| RNAseA | 200 µg/ml | 100µl from a 100mg/ml stock |

adjust the volume with MilliQ water to 50 ml

P2 buffer – can be stored at RT for up to 12 months

| Reagent | Concentrations | For 50 ml |
|---------|----------------|--------------------|
| NaOH | 200 mM | 10 ml of 1 M stock |
| SDS | 1% | 0.5 g |

adjust the volume with MilliQ water to 50 ml

P3 buffer – can be stored at RT for up to 12 months

| Reagent | Concentrations | For 50 ml |
|-------------------|----------------|-----------|
| Potassium acetate | 2.3 M | 11.3 g |

Add water till you have around 30 ml solution. Adjust the pH of the solution to **4.8** using **acetic acid**. Afterwards adjust the volume with MilliQ water to 50 ml

PB buffer – can be stored at 4 °C for up to 12 months

| Reagent | Concentrations | For 50 ml |
|-----------------------------|----------------|--|
| NaCl | 2.5 M | 25 ml of 5 M stock |
| Tris-HCl pH 8.0 | 10 mM | 0.5 ml of 1 M stock |
| EDTA | 1 mM | 0.1 ml of 0.5 M stock |
| PEG 8000 | 20%(w/v) | 10g |
| Tween 20 | 0.05% | 0.25 ml of 10% stock |
| Sera-Mag carboxylated beads | 1.14% (v/v) | 570 µl of stock (wash first, see below) |

adjust the volume with MilliQ water to 50 ml

PEG buffer – can be stored at 4 °C for up to 12 months

| Reagent | Concentrations | For 50 ml |
|-----------------|----------------|-----------------------|
| NaCl | 2.5 M | 25 ml of 5 M stock |
| Tris-HCl pH 8.0 | 10 mM | 0.5 ml of 1 M stock |
| EDTA | 1 mM | 0.1 ml of 0.5 M stock |
| PEG 8000 | 20%(w/v) | 10g |
| Tween 20 | 0.05% | 0.25 ml of 10% stock |

adjust the volume with MilliQ water to 50 ml

PE buffer – can be stored at RT for up to 12 months

| Reagent | Concentrations | For 50 ml |
|------------------------|----------------|---------------------|
| Tris-HCl pH 7.5 | 10 mM | 0.5 ml of 1 M stock |
| Ethanol | 80 % | 41.7 ml of 96% EtOH |
| MilliQ water | N/A | 7.8 ml |

EB buffer – can be stored at RT for up to 12 months

| Reagent | Concentrations | For 50 ml |
|------------------------|----------------|---------------------|
| Tris-HCl pH 8.5 | 10 mM | 0.5 ml of 1 M stock |

adjust the volume with MilliQ water to 50 ml

Consumables and equipment

Centrifuge with swing-out rotor (e.g. Rotixa 500 RS Centrifuge)

Vortex (e.g. SciQuip Vortex Mixer)







Tube rotator /shaker (e.g. Green BioResearch, Centrifuge Tube Rotator Shaker – PN: 82422201)




Magnetic stand for Eppendorf tubes (e.g. BOMB microtube racks)

Beads Washing procedure

| Step | Washing | Time | <input checked="" type="checkbox"/> |
|------|--|-------|-------------------------------------|
| 1 | Resuspend the stock of Sera-Mag carboxylated beads by vortexing. Ensure that the beads are resuspended completely before proceeding forward | 1 min | <input type="checkbox"/> |
| 2 | Take out 570 µl of beads as quickly as possible (to avoid the beads settling) and transfer them to a new microcentrifuge tube. | 1 min | <input type="checkbox"/> |
| 3 | Settle the magnetic beads on a magnetic stand and discard the supernatant. | 5 min | <input type="checkbox"/> |
| 4 | Add 1 ml of MilliQ water. Remove the tube from the magnet, resuspend the beads and centrifuge briefly. Put back the tube on the magnet and pellet the beads | 2 min | <input type="checkbox"/> |
| 5 | Remove the supernatant and add 1 ml of MilliQ. Remove the tube from the magnetic stand and resuspend the beads. After resuspension, the beads are ready to be used in preparing the PB buffer . | 5 min | <input type="checkbox"/> |

BOMB plasmid MIDI prep using Sera-Mag carboxylated beads

| Step | Clean-up | Time | <input checked="" type="checkbox"/> |
|---|--|--------|-------------------------------------|
| 1 | Pellet the bacteria from a 100 ml culture in a centrifuge at 6000 x g for 5 minutes. Discard the supernatant. | | <input type="checkbox"/> |
|  | <i>If you have a low-copy number vector, you can use up to 500 ml of culture. However, in that case you will need to multiply the volumes in the steps 2, 3, 4 by 3, as well as adjust the volume of the PEG buffer in step 8 according to the notes.</i> | | |
| 2 | Resuspend the pellet in 2 mL of P1 buffer . Make sure the pellet is completely resuspended. | | <input type="checkbox"/> |
| 3 | Add 3 ml of P2 buffer , mix by flipping the tube a few times and incubate at room temperature for 5 minutes. | 5 min | <input type="checkbox"/> |
| 4 | Add 2.5 ml of P3 buffer , mix by flipping the tube a few times. | 1 min | <input type="checkbox"/> |
|  | <i>Ensure a complete neutralization of the sample</i> | | |
| 5 | Centrifuge for 10 minutes at >14.000 x g | 10 min | <input type="checkbox"/> |
| 6 | Transfer 7 ml of the supernatant to a new tube. | 1 min | <input type="checkbox"/> |
|  | <i>Pipette carefully to avoid transferring the cell debris/precipitate. Transferring cell debris can cause protein contamination.</i> | | |
| 7 | Resuspend the beads in the PB buffer by vortexing. Make sure that all the beads are resuspended before proceeding forward. | 1 min | <input type="checkbox"/> |
|  | <i>Ensure the beads are fully resuspended otherwise you might lose your DNA.</i> | | |
| 8 | Add 3 ml of PB buffer and 4 ml of PEG buffer , mix by flipping the tube a few times. Incubate on a shaker for 30 min at RT. | 30 min | <input type="checkbox"/> |
|  | <i>Make sure that always the ratio of the recovered supernatant in step 6 and the sum of PB and PEG buffer volume is 1:1. Otherwise the DNA might not bind the beads! If you need to adjust the volumes, always adjust the PEG volume not the PB.</i> | | |
| 9 | Centrifuge the tube briefly (2-3 seconds). Settle the magnetic beads on a magnetic stand and discard the supernatant. | 5 min | <input type="checkbox"/> |
|  | <i>Ensure that the beads are completely pelleted</i> | | |
| 10 | Remove the tube from the magnetic stand and add 5 ml of PE buffer . Mix well (by pipetting, or vortex and spin down). | 2 min | <input type="checkbox"/> |
| 11 | Settle the magnetic beads on a magnetic stand and discard the supernatant. | 2 min | <input type="checkbox"/> |
| 12 | Remove the tube from the magnetic stand and add 1 ml of PE buffer . Mix well (by pipetting, or vortex and spin down). | 2 min | <input type="checkbox"/> |
| 13 | At this stage, transfer the mixture in a new microcentrifuge tube for easier work. | 1 min | <input type="checkbox"/> |
| 14 | Settle the magnetic beads on a magnetic stand and discard the supernatant. | 2 min | <input type="checkbox"/> |
| 15 | While the tube is on the magnet, add 1 ml of PE buffer . Wait for 30 seconds and then remove the supernatant. | | |
| 16 | Remove the tube from the magnetic stand and add 500 µl of PE buffer . Mix well (by pipetting, or vortex and spin down). | 2 min | <input type="checkbox"/> |
| 17 | Remove the tube from the magnet and centrifuge it briefly to pellet all the liquid. | 1 min | <input type="checkbox"/> |
| 18 | Put back the tube on the magnet and using a 10 µl pipette, remove all the remaining liquid with care to not take up any beads. | 1 min | <input type="checkbox"/> |

| | | | |
|---|--|--------|--------------------------|
| 19 | Put the tube on a heat block and dry the beads for approximately 10 minutes at 37°C. | 5 min | <input type="checkbox"/> |
|  | <i>The beads are dry when you see the beads pellet starting to crack. Do not over-dry. For midiprep, it might be that even in 10 minutes, the beads don't dry so leave them more if needed.</i> | | |
| 20 | Add 300 µl of EB buffer to elute the DNA from the beads. Pipette well to resuspend and make sure that all the beads from the side of the tube are resuspended. Incubate for 15 min at 40°C. | 15 min | <input type="checkbox"/> |
|  | <i>If you have bigger plasmid (>6kbp) it is recommended to incubate the elution at 50°C. Do not elute above 50°C.</i> | | |
| 21 | Settle the magnetic beads on a magnetic stand and transfer the supernatant to a new tube. | 5 min | <input type="checkbox"/> |
|  | <i>Make sure the beads are completely settled before transferring the supernatant. Avoid transferring any beads. If possible, leave some supernatant in the tube to minimize carryover.</i> | | |

Exemplary results




Fig 1: Comparison of commercial kits and BOMB protocols #5.3 and #5.4. MW: Gene Ruler DNA Ladder (Thermo).